## Cleaning and Sanitizing

1. Wash cleaning surfaces with solution appropriate for areas with food contact.
2. Prior to initial use, all Dr. Brown's<sup>®</sup> products should be disassembled and thoroughly cleaned with a dishwashing detergent in water and sterilized/sanitized in a microwave steam bag/electric sterilizer/top-rack dishwasher or steam sterilization (autoclave\*). Dr. Brown's<sup>®</sup> products should be washed with hot, soapy water and rinsed well between use.
3. Use Dr. Brown's<sup>®</sup> Non-Metal Bottle Brush or AccuBrush<sup>™</sup> silicone bristles to eliminate debris from the interior of the Dr. Brown's<sup>®</sup> bottle/vessel and collar.
4. Use the Wire Cleaning Brushes to thoroughly clean the internal vent system (i.e. insert and reservoir tube) and Infant-Paced Feeding Valve (blue one-way valve) only. **Hospital use is recommended to replace the wire cleaning brush every 24 hours after use.**
5. Use the chevron tip at the bottom of either bottle brush to eliminate debris inside the tip of the nipple/teat and assure the nipple/teat hole is clean and open.
6. Prior to each use and after cleaning, pull on the nipple/teat in all directions to make sure it does not tear. **Replace the nipple/teat prior to feeding if the integrity of the nipple/teat is compromised.**
7. The Dr. Brown's<sup>®</sup> bottle system, nipples/teats, brushes and one-piece silicone pacifiers may be sanitized using Dr. Brown's<sup>®</sup> Microwave Steam Sanitizing Bags.
8. Disassemble bottle system and place all cleaned parts including bottle/vessel, nipple/teat, collar, internal vent system, brushes and one-piece silicone pacifier in the bag. Place two ounces of water into the bag. **Note: Do not place wire cleaning brush into the Dr. Brown's<sup>®</sup> Microwave Steam Sanitizing Bag due to metal in the brush.**
9. Each Dr. Brown's<sup>®</sup> Microwave Steam Sanitizing bag may be used for a total of 20 cleaning cycles following the instructions for use on the bag. **Be aware of the "Warning-Safety Note" to avoid injury to hands/face.**

\*Autoclave (Steam Sterilization) is recommended not to exceed 30 cycles and is NOT recommended for the Infant-Paced Feeding Valve (blue one-way valve), Non-Metal Brush or AccuBrush<sup>™</sup>.

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