

Pacidose™ and Silicone Bulb Cleaning Brush

The Pacidose™ allows for fuss-free, accurate administration of small liquid boluses by way of a small silicone bulb with narrow tubing through the center.

- Uses include medication/vitamins, therapeutic tastings, and colostrum
- Allows infant to extract liquid in a more natural and supportive way for safe swallowing
- Attaches to a slip tip and/or catheter tip syringes
- BPA free; Individually wrapped with wire cleaning brush
- Silicone bulb is top rack dishwasher, sanitizer, and microwave steam bag compatible
- Use Silicone Bulb Cleaning Brush to clean inside channel of Pacidose™

- | | |
|--------------------------|---------------------------------|
| ✓ Mother/Baby | ✓ OR departments |
| ✓ Cardiac Units | ✓ GI/X-ray labs |
| ✓ PICU/NICU | ✓ Outpatient/Early Intervention |
| ✓ Emergency Department | ✓ Surgery Centers |
| ✓ Pediatrician's offices | ✓ NAS units |
| ✓ OB/GYN offices | ✓ Swallow Studies |

MED BULK PACKAGE

ITEM #	ITEM DESCRIPTION	RETAIL UPC	# PER CASE
HG108-MED	Pacidose™, bulb and cleaning brush only	0-72239-32255-0	25
HG107-MED	Pacidose™, bulb and cleaning brush only	0-72239-32254-0	100
SR108-MED	Silicone Bulb Cleaning Brush Only	0-72239-32621-3	240

MASTER CARTON: # UNITS VARIES/CASE

	DIMENSIONS			WEIGHT		
	WIDTH	LENGTH	HEIGHT	Cubic Ft	POUNDS	OUNCES
HG108-MED	4.0"	8.7"	6.7"	0.13 cu ft	0.0 lb	13.7 oz
HG107-MED	7.3"	12.0"	7.8"	0.40 cu ft	3.0 lb	0.5 oz
SR108-MED	3.188"	5.5"	2.125"	0.02 cu ft	0.0 lb	9.45 oz

